

SoftWick™

CoolMax® Stretch Prosthetic

HALF SOCK™

The ideal solution for managing distal volume changes with a liner!

Why pay for a whole sock when all you need is a half sock?

If you have ever received that call late on a Saturday night, you will never forget it. You cut off an old sock and put a hole in the distal end so that you could fill in some volume distally for your patient's liner. Then the sock unraveled late on a Saturday night and your patient could not get out of their prosthesis, so they called you. Sound familiar?

The Soft Wick Prosthetic Half Sock from Comfort Products is the solution to this problem with managing distal volume. The distal hole is knit right into the sock to ensure a proper fit, without unraveling, every time. The CoolMax and Lycra construction stretches to accommodate a multitude of shapes. Four sizes ensures a proper fit for almost any application. The color band knit into each sock allows easy size identification.



Soft Wick Half Sock Sizes and Part Numbers					
Size	Band Color	Length	Top Width	Toe Width	Part Number
1	Black	3"	3"	2 1/2"	SW213K
2	Gray	3"	3 1/2"	3"	SW223K
3	Royal Blue	6"	4"	2 3/4"	SW233K
4	Turquoise	6"	5 1/2"	3 1/2"	SW243K

Also available by Special Order without a hole!

Developed and manufactured by:



800-822-7500

info@comfortoandp.com